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INSTRUCTION FOR USE Please read carefully before using

PARTS DESCRIPTION	
	
SPECIFICATIONS	

Range: 32.0 °C to 43.9 °C; below 32.0 °C display L °C; above 43.9 °C display H °C
Intended use: Measurement of body temperature
Resolution: 0.1 °C
Accuracy: +/- 0.1 °C (34.0 °C to 42.0 °C), +/- 0.2 °C (other range)
Display: Liquid crystal display 3 1/2 digits
Battery (included): Micro alkaline 192 or LR41 1.55 V
Power Consumption: 0.15 milliwatt in measurement mode
More than 200 hours of continuous operation
Battery Life: IP 27(Protected against solid foreign objects with a diameter of 12.5 mm and protection against temporary immersion)
Dimensions: 121mm x 18mm x 9.8 mm
Weight: Approx. 11grams including battery
Beeper: Approx. 8 secs. sound signal when no relevant increase in temperature is detected
Memory: Built-in memory enabling display of previously measured value.
Working Conditions: Temperature: 10 °C to 40 °C; Humidity: 15 % to 95 % R.H.
 It is only possible to measure temperatures that are above the ambient temperature.

Storage Conditions: Temperature: -25 °C to +55 °C; Humidity: 15 % to 95 % R.H.
Guarantee of Quality: DIN EN ISO 80601-2-56 Medical electrical equipment. Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement (ISO 80601-2-56). The Geratherm Medical AG is certified according to DIN EN ISO 13485 and complies with Council Directive 93/42/EEC. This entitled to affix the **CE0197** marking (Notified body: TÜV Rheinland LGA Products GmbH).

Electromagnetic compatibility (EMC): The present medical product is a device with extremely high degree of noise immunity. Therefore, the enclosure with the EMC tables was not printed out.

2. PRECAUTIONS

- Use of this instrument is not intended as a substitute for consultation with your physician. Prolonged high fever requires medical attention especially for your children. Please contact your physician.
- Child's temperature shall be taken under adult's surveillance.
- Keep out of the reach of unsupervised children. Choking resulting from a child swallowing a small part that has become detached from the thermometer.
- Do not boil, bite, bend, drop or disassemble this thermometer.
- Do not attempt to disassemble the unit except to replace the battery.
- Avoid taking temperature in 30 minutes after physical exercises, bathing or dining.
- Do not store, use and transport the unit under direct sunlight, at a high temperature, or in high humidity or dust. Performance may be degraded.
- Stop using the thermometer if it operates erratically or if the display malfunctions.
- If the thermometer has been stored at below-freezing temperatures, allow the thermometer to warm naturally to room temperature before using.

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 Discard old battery carefully, out of reach of young children. Do not dispose of the battery in a fire. It may explode.

HOW TO USE

- Press the ON/OFF button. The display will show the following **1888** .
- Place the probe well under the patient's tongue. Instruct the patient to keep their mouth closed while the temperature is reading. A normal temperature by this method is usually considered between 36.0 °C and 37.5 °C.
- Do not use water!
- Place a new 1.55 VDC button size battery type 192, LR41 or equivalent in the chamber with the positive (+) side up and the negative (-) side down.
- Replace and fasten the battery cover securely.
- STOP if you meet any resistance. The sensing device is on the very tip of the probe and there is no need to insert the probe deep into the rectum. In general, the rectal temperature is around 0.5 °C to 1.0 °C **higher** than oral temperature.
- To check memory: Press on/off button and hold button for 2 seconds when the unit is under off state. The last temperature will be recalled on the LCD window.

NOTE: You may place a sanitary plastic probe shield on the sensor tip to help prevent cross-infection and contamination. Use of the probe may result in a 0.1 °C difference from actual temperature. Dispose of any used probe shield and store the unit in a safe place.

MEASUREMENT METHODS

 It is only possible to measure temperatures that are above the ambient temperature.

ORAL USE (about 45 – 90 seconds)

Place the probe well under the patient's tongue. Instruct the patient to keep their mouth closed while the temperature is reading. A normal temperature by this method is usually considered between 36.0 °C and 37.5 °C.

DO NOT USE WATER!

RECTAL USE (about 30 – 60 seconds)

Rectal temperature is recommended by your doctor, you may use a probe cover. Insert thermometer into the probe cover. Lubricate the probe cover with a water soluble jelly for easier insertion. Do not use petroleum jelly. Insert tip of probe no more than 15 mm into rectum.

STOP if you meet any resistance. The sensing device is on the very tip of the probe and there is no need to insert the probe deep into the rectum. In general, the rectal temperature is around 0.5 °C to 1.0 °C **higher** than oral temperature.

AXILLARY USE

Wipe armpit with a dry towel. Place probe in the patient's armpit and keep the patient's arm pressed firmly against their side. In general, the axillary temperature is around 0.5 °C to 1.0 °C **lower** than oral temperature. Note that in order to achieve better axillary temperature measurement result, a **minimum measuring time of 2 to 3 minutes is recommended regardless of the beep sound**.

 If you compare temperatures between two thermometers, it is possible to have two different readings even if both thermometers are within their accuracy tolerances. Differences in how you take each reading and the thermometer's tolerance range can account for nearly 0.1 °C variability.

REMEMBER: Find your normal temperature and watch for changes.

CLEANING AND DISINFECTION

- For domestic use clean the thermometer probe before and after using.
- A solution of mild detergent and cool water or 70% Isopropyl alcohol is recommended.
- Don't use hot water!
- For Disinfection please use the following disinfectant solutions. Please use one of the following disinfectants for disinfecting, following the manufacturer's instructions in each case.

Disinfectant	Dilution with water	Method
Ethanol	70-80 %	wipe the probe with a cloth moistened with the disinfectant
Isopropyl Alcohol	70-90 %	wipe the probe with a cloth moistened with the disinfectant
Glutaraldehyde	0,5-2 %	immerse the probe for 30 minutes
Aklylpolyaminomethylglycine	0,05-0,2 %	immerse the probe for 30 minutes
Benzylkonium Chloride	0,1 %	immerse the probe for 30 minutes
Benzethonium Chloride	0,1 %	immerse the probe for 30 minutes
Glucanate Chlorhexidine	0,1-0,5 %	immerse the probe for 30 minutes
Sodium Hypochlorite	0,02-0,05 %	immerse the probe for 30 minutes

 **NOTE:** Any other cleaning solution/dilution or cleaning method can cause malfunction of the thermometer.

DISPOSAL

Observe the applicable regulations when disposing of the thermometer and batteries. This thermometer must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the thermometer. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop.

SYMBOL INDEX	
	Follow instructions for use
	The device must not be disposed of with household waste.
	Store between 15 % and 95 % R.H.
	Keep dry.
	Manufacturer
	Type BF applied part
	Store between -25 °C and +55 °C
	Batch code (mm/yyyy; month/year)
	Reference number
	Model number
	Geratherm Medical AG Fahrenheitstraße 1 96331 Geratal Germany www.geratherm.com

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