NADAL® FOB QUANT Test (test cassette)

REF 272059NBUL-20

de	Gebruiksaanweisung	2	Symbols	23
en	Instructions for use	6	Our Teams	24
it	Istruzioni per l'uso	10		
pl	Sposób użycia	14		
cs	Návod k použití	18		

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16. Referenzen

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) Journal of the Canadian Association of Gastroenterology, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. ICSH; Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1986) and specifications for international haemoglobinocyanide reference preparation (3rd edition); Clin. Lab. Haemat. (1987); 9: 73-79
8. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician.(2012); 58(9): 939-948
9. Fraser C.G., et al.; Newer Fecal Tests: Opportunities for Professionals in Laboratory Medicine; Clin Chem (2012); 58(6): 963-65 DOI: 10.1373/clinchem.2012.185025
10. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65–SE87
11. Van Rossum, L.G., et al.; Cutoff value determines the performance of a semi-quantitative immunochemical faecal occult blood test in a colorectal cancer screening programme Br J Cancer (2009) 101(8) :1274-81
12. Gies A., et al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening; Gastroenterology; (2018); 154: 93-104

Rev.0, 2020-08-17 LeRe/OM

1. Intended Use

The NADAL® FOB QUANT Test is a lateral flow chromatographic immunoassay (immunological faecal occult blood test, iFOBT) for the quantitative detection of human haemoglobin in human faecal specimens. The test is intended for use as an aid in the diagnosis of pathologies of the lower gastrointestinal tract, especially in patients with abdominal complaints of unclear origin or abnormalities in bowel movements (see section 12 'Limitations'). In addition, the NADAL® FOB QUANT Test is suitable for screening asymptomatic patients who are classified as belonging to a risk group due to their age, a family history of colorectal carcinomas or other factors. The quantitative measurement and evaluation of the test result is to be carried out using the nal von minden Colibri. The test procedure is not automated and requires no special training or qualification. The NADAL® FOB QUANT Test is designed for professional use only.

2. Introduction and Clinical Significance

Colorectal cancer is one of the most commonly diagnosed types of cancer and one of the leading causes of cancer-related deaths. Nearly 55% of all colorectal cancer cases occur in more developed regions. Incidence varies worldwide by up to a factor of 10. The highest age-standardised rates can be found in Australia/New Zealand (ASR 44.8 and 32.2 per 100,000 men and women, respectively) and the lowest in West Africa (ASR 4.5 and 3.8 per 100,000 men and women, respectively).¹ Screening for occult ('hidden') blood in stool specimens can increase the chances of the early detection of colorectal cancer and thus reduce mortality rates.^{2,3}

Previous, commercially available FOB tests were based on the Guaiac method, which requires a special diet prior to testing to avoid inaccurate results caused by, for example, haem or peroxidases found in some foods. For this reason, immunological tests were developed for the specific detection of human haemoglobin in faecal specimens. Several studies have demonstrated that immunochemical methods are superior to the Guaiac method in terms of specificity, sensitivity and operability when detecting advanced adenomas and colorectal carcinomas in the lower gastrointestinal tract.^{4,6} Due to these advantages quantitative immunological FOB tests are meanwhile used in screening programmes in many countries.

3. Test Principle

The NADAL® FOB QUANT Test is a lateral flow chromatographic immunoassay (iFOBT) for the quantitative detection of human haemoglobin in human faecal specimens. Human haemoglobin is detected with the help of two specific antibodies to human haemoglobin, one of which is colour-labelled (active component of the detection reaction: approx. 0.4 µg/test) whilst the other is immobilised in the test line region (T) (active component of the detection reaction: approx. 0.7 µg/test).

Antibodies to human haemoglobin are pre-coated in the test line region (T) of the nitrocellulose membrane. During the test, the specimen extracted in the buffer reacts with antibodies to human haemoglobin which are conjugated to coloured particles and pre-coated onto the conjugate pad of the internal test strip. The mixture then migrates along the membrane by capillary action. If there is a sufficient amount of haemoglobin

in the specimen, a coloured line will develop in the test line region (T) of the membrane. The haemoglobin concentration in the specimen material correlates with the colour intensity of the test line. The line intensity is proportional to the concentration and increases from 20 ng/mL to 750 ng/mL. The concentration-dependent formation of the test line allows a rapid quantitative determination of haemoglobin in faecal specimens using the nal von minden Colibri.

The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® FOB QUANT test cassettes*
- 20 specimen collection tubes containing buffer (specimens diluent buffer)**
- 20 stool collectors
- 1 lot-specific RFID card
- 1 package insert

*containing the preservative sodium azide

**phosphate buffer contains preservatives sodium azide and ProClin™ <0.1%.

No hazard labelling is required according to Regulation (EC) N° 1272/2008 CLP. Concentrations are below the exemption thresholds.

5. Additional Materials Required

- nal von minden Colibri
- Timer, if the internal timer in the nal von minden Colibri is not to be used
- Specimen collection containers (optional, if an undiluted faecal specimen is required and the dilution step is to be carried out immediately before testing)

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use the test if the foil pouch is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection tube for each specimen obtained.

- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be disposed in accordance with local regulations.

8. Specimen Collection and Preparation

Store the specimen collection tube at room temperature prior to use. It is recommended to collect faecal specimen using a stool collector. Avoid dilution of faecal specimens with urine or water from the toilet bowl.


1. Hold the specimen collection tube upright and open the light-blue cap with the applicator stick (spiral rod). Remove it. Be careful not to spill or splash buffer solution from the collection tube.
2. Collect specimens by inserting the applicator stick in 3 different sites of the faecal specimen.

Note: Please be sure to insert the applicator stick three times in a row into the faecal specimen. Do NOT put the applicator stick back into the tube in between collecting the 3 specimens. Be careful not to spill fluid from the tube. The quality of the specimen and compliance with these instructions will affect your test result.

3. Place the applicator stick, along with the specimen, back into the specimen collection tube and close it properly.
4. Shake the specimen collection tube vigorously so that the specimen is fully mixed with the diluent buffer.

Note:

If a specimen is collected in the specimen collection container by the patient themselves, it is necessary to inform the patient how to collect, store and transport the faecal specimen as well as which conditions must be avoided.

The NADAL® FOB QUANT Test is intended for use with human faecal specimens only. Patients should not collect faecal

specimens during their menstrual period or 3 days before and after it, or if they have bleeding hemorrhoids, blood in the urine or strongly forced bowel movements.

Aspirin and other medications taken in excess may cause gastrointestinal irritation, resulting in occult bleeding. If possible, such substances should be discontinued - following consultation with the doctor - at least 48 hours prior to testing.

No dietary restrictions are necessary before testing. Alcohol is suspected of irritating the intestines and should be avoided for at least 48 hours prior to testing.

Urine and the excessive dilution of specimens with toilet water can lead to inaccurate test results.

Perform the test as soon as possible after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. On the whole, the storage time of specimens at 2-8°C should not exceed 7 days. Specimens can be transported to the doctor's office at room temperature without any problems; however this should be done within 2 hours. Bring specimens to room temperature prior to testing.

9. Test Procedure

I. Preparation


1. Bring foil pouches and patient specimens (extracted specimens) to room temperature (18-30°C) prior opening the pouches.
2. Remove the NADAL® FOB QUANT test cassette from its sealed pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch.
3. Switch on the nal von minden Colibri and make sure that the test-specific RFID card is at hand.

Note: Every lot has its own calibration data. This is saved on a RFID card delivered with every pack of tests. In order to correctly interpret the test, the lot-specific RFID card must be placed on the nal von minden Colibri. To avoid confusion between cards, old RFID cards should be disposed of once tests from the corresponding lot have been used up.

II. Test Procedure


1. Shake the specimen collection tube carefully to ensure that the faecal specimen is completely mixed with the dilution buffer.
2. Unscrew the rounded white cap on the collection tube. Using a piece of tissue, break off the tip of the tube with a twisting motion.

3. Holding the specimen collection tube vertically, discard the first drop and then transfer 3 drops of the specimen solution into the round specimen well (S) on the test cassette.

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

4. Start the timer or use the integrated timer in the nal von minden Colibri.
5. Place the nal von minden Colibri onto the test cassette, so that it fits exactly into the cavity.
6. Read the test result after exactly 5 minutes.

10. Result Interpretation

Following an analysis, the nal von minden Colibri displays the quantitative result on-screen.

The measuring range is 20 to 750 ng/mL (this corresponds to approx. 2 to 75 µg haemoglobin/g stool). The result is shown in the units "approx. µg/g" and in "ng/mL" on the Colibri display. Concentrations below this level are displayed as "approx. <2 µg/g" and "<20 ng/mL", whilst concentrations above this level are displayed as "approx. >75 µg/g" and ">750 ng/mL". Please select the unit that corresponds to your laboratory requirements and ensure that the concentration value on the Colibri display is shown directly before the unit.

Invalid:

The control line (C) fails to appear. In this case, the nal von minden Colibri shows "FAIL" on the display and a low-pitched acoustic signal can be heard. Results from any test in which the nal von minden Colibri has displayed "FAIL" must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Calibration

The NADAL® FOB QUANT Test was calibrated using the ClinChek® haemoglobin control (RECIPE CHEMICALS + INSTRUMENTS GmbH), the concentration of which was determined using the cyanmethemoglobin method (recommended reference method of the ICSH⁷). The calibration data is saved on the RFID chip included. It is therefore not necessary for the user to calibrate the device.

12. Quality Control

The internal procedural control is included in the NADAL® FOB QUANT Test:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Control standards are not supplied with the NADAL® FOB QUANT Test. Nevertheless, it is recommended that positive and negative controls are tested as part of *good laboratory*

practice (GLP) with every new lot or delivery and once a week or according to the laboratory standard procedure in order to confirm the test procedure and to verify proper test performance. Control materials should be tested in the same way as patient specimens. Furthermore, you can verify the correct functionality of the device using the included QC cassette. For detailed information, please consult the device manual.

13. Limitations

- The NADAL® FOB QUANT Test is for professional *in-vitro* diagnostic use only. It should be used for the quantitative detection of human haemoglobin in faecal specimens only.
- Other than colorectal bleeding, the presence of blood in the stool may have a number of different causes, such as haemorrhoids, blood in the urine or stomach irritation. Menstrual blood, minor injuries (e.g. from hard stools, rectal medication or digital rectal examinations) or diseases (e.g. diarrhoeal disease, inflammatory bowel disease) may also lead to blood in the stool. Drugs that inhibit blood clotting or irritate the intestinal mucosa are also suspected of exacerbating bleeding tendencies. However, the influence of drugs is controversial, as studies often do not demonstrate any significant differences on the positive predictive value. For a comprehensive review, see references.⁸
- Negative results do not exclude bleeding, as blood may be unevenly distributed in stool specimens. In addition, some polyps and colorectal carcinomas can bleed intermittently or not at all. In particular, colorectal polyps in the very early stages often go undetected. Some cancer precursors are detected more reliably using colonoscopy than iFOBT, which is why colonoscopy continues to be regarded as the 'gold standard'. However, as an invasive examination procedure, it poses a greater risk for the patient.
- The NADAL® FOB QUANT Test is not suitable for the detection of bleeding in the upper gastrointestinal tract, as haemoglobin degrades as it passes through the gastrointestinal passage.
- Not all colorectal bleeding is due to precancerous or cancerous polyps. As with all diagnostic tests, a confirmed diagnosis should only be made by the physician after all clinical and laboratory findings have been evaluated. Positive test results must be clarified by means of follow-up examinations. Persistent symptoms despite a negative result should be clarified by differential diagnosis.

14. Performance Characteristics

Clinical performance

In a comparison study between the NADAL® FOB QUANT Test and the RIDASCREEN® Haemoglobin ELISA, the examination of patient specimens showed a correlation of 91.3%.

It is recommended to adjust the cut-off according to the intended use (objective of the examination, e.g. screening or supportive diagnostics, national guidelines, patient age and sex).^{9,10,11} The smaller the cut-off is, the higher the sensitivity and the lower the specificity.¹²

Analytical performance

Detection limit, measuring range and hook effect:

The detection limit of the NADAL® FOB QUANT Test for the quantitative evaluation is 20 ng/mL or approx. 2 µg/g stool. The measuring range is between 20 ng/mL and 750 ng/mL,

which corresponds to approx. 2 µg/g or 75 µg/g stool and can be quantitatively evaluated using the nal of minden Colibri. No quantification is possible below 20 ng/mL and above 750 ng/mL.

No hook effect relevant to the results can be detected up to a concentration of 10,000 ng/mL.

Analytical specificity

The NADAL® FOB QUANT Test is specific for human haemoglobin and shows no cross-reaction with haemoglobin from bovine, pig, horse, sheep, goat, chicken, turkey, rabbit or fish blood at concentrations of up to 0.5 mg/mL.

The test shows no cross-reactivity to haemoglobin S. Haemoglobin F with a concentration of up to 3.4% has no effect on the test result.

Interfering substances

The NADAL® FOB QUANT Test shows no interference with the listed substances up to the specified concentrations:

Analyte	Concentration
Acetylsalicylic acid	100 mg/L
Ascorbic acid	200 mg/L
Bilirubin	50 mg/L
Caffeine	400 mg/L
Iron(II) sulfate heptahydrate	25 mg/mL
Glucose	200 mg/L
Uric acid	50 mg/L
Urea	200 mg/L
Human albumin	25 mg/L
Lipids*	20 mg/L

* Fatty acids saturated: 14.1%, monounsaturated: 76.1%, polyunsaturated: 9.8%

Precision

The repeatability and reproducibility of the test were established by testing haemoglobin concentrations of 20 ng/mL, 150 ng/mL and 750 ng/mL in a 10-fold determination with 3 independent lots. For repeatability, the mean intra-lot variance is 7.5% in the low, 4.6% in the medium and 13.9% in the high concentration range. The mean inter-lot variance for the reproducibility is 7.3% in the low, 5.0% in medium and 9.8% in the high concentration range.

15. Serious incident reporting

In the event of any serious incidents related to the performance of the NADAL® FOB QUANT Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

16. References

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.J., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) *Journal of the Canadian Association of Gastroenterology*, gwz009, <https://doi.org/10.1093/jcag/gwz009>

7. ICSH; Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1986) and specifications for international haemoglobinocyanide reference preparation (3rd edition); Clin. Lab. Haemat. (1987); 9: 73-79
8. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician. (2012); 58(9): 939-948
9. Fraser C.G., et al.; Newer Fecal Tests: Opportunities for Professionals in Laboratory Medicine; Clin Chem (2012); 58(6): 963-65 DOI: 10.1373/clinchem.2012.185025
10. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65-SE87
11. Van Rossum, L.G., et al.; Cutoff value determines the performance of a semi-quantitative immunochemical faecal occult blood test in a colorectal cancer screening programme Br J Cancer (2009) 101(8): 1274-81
12. Gies A., et al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening: Gastroenterology; (2018); 154: 93-104

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16. Bibliografia

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86

2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) *Journal of the Canadian Association of Gastroenterology*, gwz009, <https://doi.org/10.1093/cag/gwz009>
7. ICSH; Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1986) and specifications for international haemoglobinocyanide reference preparation (3rd edition); Clin. Lab. Haemat. (1987); 9: 73-79
8. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician (2012); 58(9): 939-948
9. Fraser C.G., et al.; Newer Fecal Tests: Opportunities for Professionals in Laboratory Medicine; Clin Chem (2012); 58(6): 963-65 DOI: 10.1373/clinchem.2012.185025
10. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65-SE87
11. Van Rossum, L.G., et al.; Cutoff value determines the performance of a semi-quantitative immunochemical faecal occult blood test in a colorectal cancer screening programme Br J Cancer (2009) 101(8): 1274-81
12. Gies A., et al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening; Gastroenterology; (2018); 154: 93-104

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16. Bibliografia

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal-occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal-occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.J., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) Journal of the Canadian Association of Gastroenterology, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. ICSH; Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1986) and specifications for international haemoglobin cyanide reference preparation (3rd edition); Clin. Lab. Haemat. (1987); 9: 73-79
8. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician. (2012); 58(9): 939-948
9. Fraser C.G., et al.; Newer Fecal Tests: Opportunities for Professionals in Laboratory Medicine; Clin Chem (2012); 58(6): 963-65 DOI: 10.1373/clinchem.2012.185025
10. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65-SE87
11. Van Rossum, L.G., et al.; Cutoff value determines the performance of a semi-quantitative immunochemical faecal occult blood test in a colorectal cancer screening programme Br J Cancer (2009) 101(8) :1274-81
12. Gies A., et al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening; Gastroenterology; (2018); 154: 93-104

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16. Literatura

1. Ferlay, J., et al.; Cancer incidence and mortality worldwide: Sources, methods and major patterns in GLOBOCAN 2012; Int J Cancer (2015); 136(5): E359-86
2. Hardcastle J.D., et al.; Randomised controlled trial of faecal occult-blood screening for colorectal cancer; Lancet (1996); 348 (9040): 1467-71
3. Kronborg O., et al.; Randomised study of screening for colorectal cancer with faecal occult-blood test; Lancet (1996); 348 (9040): 1467-71
4. Allison J.E., et al.; Population screening for colorectal cancer means getting FIT: The past, present, and future of colorectal cancer screening using the fecal immunochemical test for hemoglobin (FIT); Gut Liver (2014); 8 (2): 117-30
5. Park D.I., et al.; Comparison of guaiac-based and quantitative immunochemical fecal occult blood testing in a population at average risk undergoing colorectal cancer screening; Am J Gastroenterol (2010); 105 (9): 2017-25
6. Sultanian R., et al.; The impact of transitioning from guaiac-fecal occult blood testing to fecal immunochemical testing in a Canadian colon cancer screening program (2019) Journal of the Canadian Association of Gastroenterology, gwz009, <https://doi.org/10.1093/jcag/gwz009>
7. ICSH; Recommendations for reference method for haemoglobinometry in human blood (ICSH standard 1986) and specifications for international haemoglobin cyanide reference preparation (3rd edition); Clin. Lab. Haemat. (1987); 9: 73-79
8. Konrad G.; Are medication restrictions before FOBT necessary?; Can Fam Physician. (2012); 58(9): 939-948
9. Fraser C.G., et al.; Newer Fecal Tests: Opportunities for Professionals in Laboratory Medicine; Clin Chem (2012); 58(6): 963-65 DOI: 10.1373/clinchem.2012.185025
10. Halloran S.P., et al.; European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition DOI <http://dx.doi.org/10.1055/s-0032-1309791> Endoscopy 2012; 44: SE65–SE87

11. Van Rossum, L.G., et al.; Cutoff value determines the performance of a semi-quantitative immunochemical faecal occult blood test in a colorectal cancer screening programme Br J Cancer (2009) 101(8) :1274-81
12. Gies A., et. al.; Direct Comparison of Diagnostic Performance of 9 Quantitative Fecal Immunochemical Tests for Colorectal Cancer Screening; Gastroenterology; (2018); 154: 93-104

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