# NADAL<sup>®</sup> Vitamin D Quant Test (test cassette)

REF 2260004NBUL-25

<b>de</b> Gebrauchsanweisung	2	<b>cs</b> Návod k použití	25
<b>en</b> Instructions for use	6	Symbols	31
<b>fr</b> Instructions d'utilisation	9	Our Teams	32
<b>es</b> Instrucciones de uso	13		
<b>it</b> Istruzioni per l'uso	17		
<b>pl</b> Sposób użycia	21		

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**13. Referenzen**

- Holick, MF. Vitamin D status: Measurement, Interpretation and clinical application. Ann. Epidemiol. 2009, 19(2):73-78.
- Morris HA. Vitamin D: A Hormone for All Seasons – How much is enough? Clin. Biochem. Rev., 2005, 26:21-32.
- Moyad MA. Vitamin D: a rapid review. Dermatol Nurs. 2009, 21:25-30.
- Zerwekh JE. Blood biomarkers of vitamin D status. Am J. Clin Nutr. 2008, 87:1087S-915
- Schöttker B, et al. Vitamin D and mortality: meta-analysis of individual participant data from a large consortium of cohort studies from Europe and the United States. BMJ. 2014, 348:g3656
- Jenkinson C et al. High throughput LC-MS/MS method for the simultaneous analysis of multiple vitamin D analytes in serum. Chromatogr B Analyt Technol Biomed Life Sci 2016 Mar 1; 1014:56-63.

Rev.0, 2022-05-24 DH

### 1. Intended Use

The NADAL® Vitamin D Quant Test is a lateral flow test for the quantitative determination of total 25-hydroxy Vitamin D (25-OH Vitamin D) in human fingerstick blood or serum. This assay is intended for use in screening for a Vitamin D deficiency. Liquid chromatography with tandem mass spectrometry (LC-MS/MS) assays or other quantitative immunoassays are recommended to further confirm the diagnostic test results. The quantitative measurement and evaluation of the test result have to be carried out using the Colibri point-of-care reader. The test is for professional use only.

### 2. Introduction and Clinical Significance

Vitamin D is a steroid hormone responsible for enhancing the intestinal absorption of calcium and the regulation of its homeostasis. The two common forms of Vitamin D are Vitamin D2 and Vitamin D3. Vitamin D3 is naturally produced in the human skin through exposure to ultraviolet light whereas Vitamin D2 is mainly obtained from food. Vitamin D is transported to the liver where it is metabolised to 25-hydroxy Vitamin D. The 25-hydroxy Vitamin D blood test is used to determine Vitamin D concentration in the body. The blood concentration of 25-hydroxy Vitamin D is considered the best indicator of Vitamin D status.

Vitamin D deficiency is now recognised as a global epidemic. Virtually every cell in the human body has receptors for Vitamin D, meaning that they all require a sufficient Level of Vitamin D to function adequately. The health risks associated with Vitamin D deficiency are far more severe than previously thought. Vitamin D deficiency has been linked to various serious diseases: Osteoporosis, osteomalacia, multiple sclerosis, cardiovascular diseases, pregnancy complications, diabetes, depression, strokes, autoimmune diseases, flu, COVID-19, infectious diseases, various cancers, alzheimer's, obesity and higher mortality rates etc. As a result, measuring (25-OH) Vitamin D levels is now considered medically necessary, and maintaining sufficient levels is considered not just vital for improving bone health, but also overall health and well-being.

Multiple guidelines for Vitamin D deficiency have been published by various health organisations; but a universal recommendation remains to be established. Recent literature has suggested the following ranges for the classification of Vitamin D status:

25-OH Vitamin D Level	Reference range (ng/mL)	Reference range (nmol/l)
Deficiency	0-10	0-25
Insufficient	10-30	25-75
Sufficient	30-100	75-250
Excessive, but not toxic	100-150	250-375
Toxic	>150	>375

### 3. Test Principle

The NADAL® Vitamin D Quant Test allows to measure 25-OH Vitamin D quantitatively via interpretation of the test line's colour intensity. The NADAL® Vitamin D Quant Test is an immunochromatographic test based on two specific

monoclonal antibodies against 25-OH Vitamin D, one conjugated with colloidal gold and the other immobilised on the test membrane. The concentration-dependent formation of the test line coloration allows a rapid, quantitative determination of 25-OH Vitamin D in human fingerstick blood or serum.

As the sample material migrates by capillary action along the membrane within the test cassette, the coloured anti-25-OH Vitamin D-colloidal gold conjugate complexes with 25-OH Vitamin D from the sample. This complex migrates further along membrane to the test region (T), where it is captured by another anti-25-OH Vitamin D immobilized on the membrane, leading to the formation of a pink / purple coloured line. The intensity of the coloured line in the test line region correlates with the 25-OH Vitamin D concentration: the higher the concentration of 25-OH Vitamin D in the sample tested, the stronger is the colour intensity of the line. A control line is present in the test window to work as procedural control. This coloured line should always appear on the control line region (C) if the test cassette is stored under recommended storage conditions and the test is performed appropriately.

### 4. Reagents and Materials Supplied

- 25 NADAL® Vitamin D Quant Test
- 26 collection tubes containing buffer
- 26 blood collectors
- 1 RFID card
- 1 package insert

### 5. Additional Materials Required

- Lancets
- Alcohol swabs
- Colibri point-of-care reader
- Timer (unless the integrated timer in the Colibri is used)
- Vitamin D control
- 5 µl micropipette (for use with serum samples or Vitamin D control)

### 6. Storage and Stability

Test kits should be stored at 4-30°C and are stable until the expiration date stated on the foil pouches. Test cassettes must remain in sealed foil pouches until use as they are sensitive to humidity. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

### 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- The correct RFID card must be used in order to obtain reliable results. RFID cards are specific to tests and lots.
- Do not add specimens to the reaction area (result area).

- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new collection tube and a new blood collector for each specimen.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

## 8. Specimen Collection and Preparation

### Preparation

Before performing the test, please make sure that all components are brought to room temperature (see "Storage and Stability").

### Obtaining a blood sample

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.
- Collect 10 µl blood using the blood collector.

## 9. Test Procedure

### Testing a fingerstick blood sample

**Bring tests, samples, buffer and/or controls to room temperature prior to testing.**


1. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
2. Place the test cassette on a clean and level surface.
3. Tap the collection tube on a flat surface to ensure that all of the liquid is at the bottom of the tube.

4. Peel off the aluminium foil from the top of the collection tube.
5. Collect blood according to the instructions above ("Obtaining a blood sample"). Gently touch the tip of blood collector to the blood droplet. Capillary action will draw up 10 µL of blood into the blood collector and stop automatically.
6. Check the level of blood in the blood collector. The capillary inside the collector must be completely filled with blood.
7. Fully insert the blood collector into the collection tube and push firmly to close tightly.
8. Shake the combined tube with a thrusting motion 3-4 times to transfer the blood from the blood collector into the collection tube and shake it until the solutions are mixed completely.
9. Remove the cap of the blood collector.
10. Invert the tube with the blood sample and gently squeeze it to apply 3 drops into the sample well (S) of the test cassette.
11. Read and record the results after 15 minutes using the Colibri point-of-care reader (for use of the Colibri point-of-care reader please refer to the corresponding manual).

### Important Notes:

- Temperature can adversely affected test results. Tests should therefore be carried out at room temperature (21-24°C if possible).
- For optimum results, ensure the fingerstick blood sample is completely mixed with the buffer. This can be determined by checking the uniformity of the red colouring of the pre-mixed blood in the collection tube and blood collector. Insufficient mixing can cause incorrect test results.
- Once mixed, the collection tube should be squeezed **gently** in order to dispense three full sample drops into the sample well (S).
- Results should be interpreted using the Colibri point-of-care reader **exactly 15 minutes** after dispensing the sample into the sample well (S). Failure to do so may lead to inaccurate results.

- Anticoagulated blood or plasma samples must not be used for testing with the NADAL® Vitamin D Quant Test, as anticoagulants will influence the test results.

#### Testing with vitamin D control / serum:

The NADAL® Vitamin D Quant Test for use with fingerstick blood. However, the vitamin D control or serum samples can also be used. For test preparation, refer to step 1-4 of 'Testing a fingerstick blood sample' (see above). Instead of obtaining a fingerstick sample with the blood collector, apply 5 µL of Vitamin D control or serum into the collection tube using a micropipette (not provided with the test kit) and proceed with step 7-11 of 'Testing a fingerstick blood sample' (see above).

#### Important note:

- Make sure exactly 5 µL serum/control solution is transferred into the collection tube. If there is liquid on the outside of the pipette tip, carefully wipe it off with a paper towel in the direction of the pipette tip. Due to the small sample volume, even the smallest changes in volume can lead to an inaccurate measurement result.

#### Result interpretation:

The Colibri point-of-care reader automatically calculates the final result based on the detected intensity of the test line and the calibration curve stored on the corresponding RFID card. 25-OH Vitamin D levels are measured in ng/mL. Results obtained with the NADAL® Vitamin D Quant Test are interpreted using the following categories:

25-OH Vitamin D Level (ng/mL)	Categorisation
< 10.0 ng/mL	DEFICIENCY
≥ 10.0 – < 30.0 ng/mL	INSUFFICIENCY
≥ 30.0 – < 100.0 ng/mL	SUFFICIENT
≥ 100.0 ng/mL	EXCESSIVE OR TOXIC

#### 10. Quality Control

An internal procedural control is included in the test cassette: A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking. Good laboratory practice recommends the use of external control materials to ensure proper test kit performance.

If control values do not fall within the established range, assay results should be considered invalid. A set of two LC-MS/MS-confirmed vitamin D controls can be purchased separately.

#### 11. Performance Characteristics

##### Sensitivity:

The sensitivity of the NADAL® Vitamin D Quant Test is 3 ng/mL (7.5 nmol/L). The sensitivity was determined by calculating the mean plus 3.3 times of standard deviation of 20 vitamin D-free serum tests.

##### Detection range:

The detection range of the NADAL® Vitamin D Quant Test with the Colibri point-of-care reader is from 3 ng/mL (7.5 nmol/L) to 100 ng/mL (250 nmol/L).

##### Accuracy:

The accuracy of the NADAL® Vitamin D Quant Test was evaluated using human fingerstick blood samples in comparison with a reference 25OH Vitamin D ELISA assay using corresponding serum samples. The comparative result

showed a linear regression with a slope of 1.02 and a correlation coefficient of 92%. In conclusion, the NADAL® test results using human fingerstick blood samples showed good agreement with the ELISA results using corresponding serum samples. The accuracy of the NADAL® Vitamin D Quant Test was also evaluated using 20 serum samples in comparison with LC-MS/MS Assay (the standard for 25-OH Vitamin D measurement). The comparison result showed a linear regression with a slope of 0.98 and a correlation coefficient of 98%. In conclusion, the NADAL® Vitamin D Quant Test results agree closely with the true values generated using an LC-MS/MS assay.

#### Precision:

##### Intra Lot

Sample	No. of Replicates	Mean (ng/mL)	Standard Deviation (ng/mL)	Coefficient Variation (%)
Serum	20	40.4	3.1	7.6%
Blood-1	10	33.0	2.5	7.6%
Blood-2	10	38.7	3.9	10.1%

##### Inter Lot

Sample	No. of Lot	No. of Replicates	Mean (ng/mL)	Standard Deviation (ng/mL)	Coefficient Variation (%)
Serum #1	3	60	33.0	4.7	14.3%
Serum #2	3	60	55.9	6.1	11.0%
Blood-1	3	15	38.8	4.3	11.8%

#### Specificity:

30 serum samples free from vitamin D were tested and all showed negative results; indicating 100% specificity.

No interference or cross reactivity was observed with bilirubin, triglycerides, cholesterol, vitamin B12 and vitamin C.

#### 12. Expected Results

The NADAL® Vitamin D Quant Test is a rapid quantitative assay. The test is intended for use in screening individuals to ascertain vitamin D levels. This assay only provides a preliminary analytical test result. The liquid chromatography with tandem mass spectrometry (LC-MS/MS) assays or quantitative immunoassays are recommended to confirm the analytical result.

#### 13. References

1. Holick, MF. Vitamin D statuses: Measurement, Interpretation and clinical application. Ann. Epidemiol. 2009, 19(2):73-78.
2. Morris HA. Vitamin D: A Hormone for All Seasons – How much is enough? Clin. Biochem. Rev., 2005, 26:21-32.
3. Moyad MA. Vitamin D: a rapid review. Dermatol Nurs. 2009, 21:25-30
4. Zerwekh JE. Blood biomarkers of vitamin D status. Am J. Clin Nutr. 2008, 87:10875-915
5. Schöttker B, et al. Vitamin D and mortality: meta-analysis of individual participant data from a large consortium of cohort studies from Europe and the United States. BMJ. 2014, 348:g3656
6. Jenkinson C et al. High throughput LC-MS/MS method for the simultaneous analysis of multiple vitamin D analytes in serum. Chromatogr B Analyt Technol Biomed Life Sci 2016 Mar 1; 1014:56-63.

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**13. Referencias**

- Holick, MF. Vitamin D status: Measurement, Interpretation and clinical application. Ann Epidemiol. 2009, 19(2):73-78.
- Morris HA. Vitamin D: A Hormone for All Seasons – How much is enough? Clin. Biochem. Rev. 2005, 26:21-32.
- Moyad MA. Vitamin D: a rapid review. Dermatol Nurs. 2009, 21:25-30
- Zerwekh JE. Blood biomarkers of vitamin D status. Am J. Clin Nutr. 2008, 87:1087S-91S
- Schöttker B, et al. Vitamin D and mortality: meta-analysis of individual participant data from a large consortium of cohort studies from Europe and the United States. BMJ. 2014, 348:g3656
- Jenkinson C et al. High throughput LC-MS/MS method for the simultaneous analysis of multiple vitamin D analytes in serum. Chromatogr B Analyt Technol Biomed Life Sci 2016 Mar 1; 1014:56-63.

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**13. Bibliografia**

- Holick, MF. Vitamin D status: Measurement, Interpretation and clinical application. Ann. Epidemiol. 2009, 19(2):73-78.
- Morris HA. Vitamin D: A Hormone for All Seasons – How much is enough? Clin. Biochem. Rev., 2005, 26:21-32.
- Moyad MA. Vitamin D: a rapid review. Dermatol Nurs. 2009, 21:25-30
- Zerwekh JE. Blood biomarkers of vitamin D status. Am J. Clin Nutr. 2008, 87:1087S-91S
- Schöttker B, et al. Vitamin D and mortality: meta-analysis of individual participant data from a large consortium of cohort studies from Europe and the United States. BMJ. 2014, 348:g3656
- Jenkinson C et al. High throughput LC-MS/MS method for the simultaneous analysis of multiple vitamin D analytes in serum. Chromatogr B Analyt Technol Biomed Life Sci 2016 Mar 1; 1014:56-63.

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#### 13. Reference

- Holick, MF. Vitamin D statuses: Measurement, Interpretation and clinical application. Ann. Epidemiol. 2009, 19(2):73-78.
- Morris HA. Vitamin D: A Hormone for All Seasons – How much is enough? Clin. Biochem. Rev. 2005, 26:21-32.
- Moyad MA. Vitamin D: a rapid review. Dermatol Nurs. 2009, 21:25-30
- Zerwekh JE. Blood biomarkers of vitamin D status. Am J. Clin Nutr. 2008, 87:10875-915
- Schöttker B, et al. Vitamin D and mortality: meta-analysis of individual participant data from a large consortium of cohort studies from Europe and the United States. BMJ. 2014, 348:g3656
- Jenkinson C et al. High throughput LC-MS/MS method for the simultaneous analysis of multiple vitamin D analytes in serum. Chromatogr B Analyt Technol Biomed Life Sci 2016 Mar 1; 1014:56-63.


Rev.0, 2022-05-24 JV

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico- diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur- begränsning	Temperatuurlimiet	Temperatur- begrænsning	Temperatur begrænsning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merkning
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for<n> tester

## Our Teams

### Germany:

#### Regensburg

Tel: +49 941 290 10-0  
Fax: +49 941 290 10-50

#### Moers

Tel: +49 2841 99820-0  
Fax: +49 2841 99820-1

#### Austria:

Tel: +49 941 290 10-29  
Free Tel: 0800 291 565  
Fax: +49 290 10-50  
Free Fax: 0800 298 197

#### UK & Ireland:

Tel: +49 941 290 10-18  
Free Tel –UK: 0808 234 1237  
Free Tel – IRE: 1800 555 080  
Fax: +49 290 10-50

#### France:

France Tel: 0800 915 240  
France Fax: 0800 909 493

#### Switzerland

Swiss Tel: 0800 564 720  
Swiss Fax: 0800 837 476

#### Belgium

Belgium Tel: 0800 718 82  
Belgium Fax: 0800 747 07

### Luxembourg

Lux, Tel: 800 211 16  
Lux, Fax: 800 261 79

### Spain:

Tel: +49 941 290 10-759  
Free Tel: 900 938 315  
Fax: +49 941 290 10-50  
Free Fax: 900 984 992

### Italy:

Tel: +49 941 290 10-34  
Fax: +49 941 290 10-50

### Poland:

Tel: +49 941 290 10-44  
Free Tel: 00 800 491 15 95  
Fax: +49 941 290 10-50  
Free Fax: 00 800 491 15 94

### Portugal:

Tel: +49 941 290 10-735  
Tel, Verde: 800 849 230  
Fax: +49 941 290 10-50  
Fax Verde: 800 849 229

### Netherlands:

Tel: +31 30 75 600  
Free Tel: 0800 0222 890  
Fax: +31 70 30 30 775  
Free Fax: 0800 024 9519

### Nordic countries:

#### Denmark

Tel: +31 703075 605  
Free Tel: 808 887 53

#### Finland

Tel: +31 703075 606  
Free Tel: 0800 918 263  
Free Fax: 0800 918 262

#### Norway

Tel: +31 703075 605  
Free Tel: 800 16 731

#### Sweden

Tel: +31 703075 605  
Free Tel: 020 79 09 06

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Tel: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH