# NADAL<sup>®</sup> Influenza A/B Scan Test (test cassette)

REF 242006NBUL-10

<b>de</b>	Gebrauchsanweisung	2	Symbols	16
<b>en</b>	Instructions for use	8	Our Teams	16

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### 15. Referenzen

- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

Rev. 4, 2021-01-20 KBe/SDe

#### 1. Intended Use

The NADAL® Influenza A/B Scan Test is a rapid chromatographic immunoassay for the qualitative detection of influenza virus type A and B antigens (nucleoproteins) extracted from nasal/nasopharyngeal swabs and nasal washes/aspirates of symptomatic patients. The test is intended as an aid in the rapid, differential diagnosis of influenza virus A and B infections. It is not intended for the detection of influenza C viral antigens. The qualitative measurement and evaluation of the test result is to be carried out using the nal von minden Colibri Point of Care reader. Negative results do not exclude influenza virus A or B infection and should be confirmed by cell culture or a molecular biological assay. The test is designed for professional use only.

#### 2. Introduction and Clinical Significance

Influenza is a highly contagious viral infection of the upper respiratory tract which is characterised by antigen variability, seasonality and the impact on the general population. Of the two main types (A and B) of influenza viruses, Influenza A subtypes are differentiated by antigen variability of the surface glycoproteins (hemagglutinin and neuraminidase). Influenza A virus is the most prevalent and is associated with the most serious epidemics. Influenza can cause severe complications such as bronchitis or pneumonia, particularly in children, elderly people or people with chronic respiratory disease. Yet, it most commonly occurs as a mild viral infection transmitted by respiratory secretions through sneezing or coughing. There are many other viral infections that can mimic influenza symptoms, making laboratory tests necessary to distinguish it from other acute respiratory infections. New effective antivirals are available since the late 1990s but these drugs are effective if administered early (before 48 hours after the onset of the illness). Virus isolation is still considered as the gold standard method for influenza diagnosis, with a sensitivity of nearly 100% after 3 days. However, patient health care and economic costs could be greatly improved by the use of rapid, specific and sensitive tests for antigen detection in order to enable a successful antiviral treatment.

#### 3. Test Principle

The NADAL® Influenza A/B Scan Test enables the detection of influenza virus type A and B antigens through analysis of colour development on the test strip with the help of the nal von minden Colibri Point of Care reader. Anti-influenza type A and B antibodies to nucleoprotein antigens are immobilised in the test line region of the membrane respectively. During the test, extracted antigens bind to further anti-influenza type A and B antibodies conjugated to coloured particles and pre-coated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there are sufficient influenza virus type A or B antigens in the specimen, a coloured line will form in the test line region 'A' or 'B' of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result.

The appearance of a coloured line in the control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

#### 4. Reagents and Materials Supplied

- 10 NADAL® Influenza A/B Scan test cassettes
- 10 extraction tubes incl. dropper-caps
- Additional material provided according to 93/42/EEC:  
Due to possible supply shortages of accessory medical products as a result of the current COVID-19 pandemic, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.

a) 10 sterile swabs CE 2797

Puritan Medical Products Company LLC  
31 School Street  
Guilford, Maine 04443-0149 USA (authorised EU representative EMERGO EUROPE, The Hague, The Netherlands)

b) 10 sterile swabs, CE 0197

Jiangsu Changfeng Medical Industry Co., Ltd  
Touqiao Town, Guangling District, Yangzhou,  
Jiangsu 225109 China (authorised EU representative Lins Service & Consulting GmbH, Obere Seegasse 34/2, 69124 Heidelberg, Germany)

- 1 Buffer (7 mL), contains sodium azide (<0.1%) as preservative
- 1 reagent holder
- 1 package insert

#### 5. Additional Materials Required

- nal von minden Colibri Point of Care reader
- Stopwatch if internal timer of nal von minden Colibri is not used
- If necessary: physiological saline solution, syringes and specimen containers (for nasal washes)
- If necessary: an aspirating tube with a trap (for nasal aspirates)
- If necessary: 300 µL pipettes (for nasal washes/aspirates)

#### 6. Storage & Stability

The NADAL® Influenza A/B Scan test kits should be stored at 2-30°C and used before the expiry date printed on the packaging. The test cassette must remain in the sealed foil pouch until use. Do not freeze the test.

#### 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure of test and used devices prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits. Do not mix dropper-caps.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.

- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Do not use swabs from damaged pouches.
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

## 8. Specimen Collection and Preparation

### Specimen Collection

Acceptable specimens for testing with the NADAL® Influenza A/B Scan Test include nasal/nasopharyngeal swabs as well as nasal washes/aspirates. Do not use specimens that are obviously contaminated with blood as it may interfere with the flow of sample on the test strip and lead to inaccurate test results. Use freshly collected specimens for the best test performance. Rapid tests have more reliable clinical performance when performed early in the course of infection.<sup>4</sup>

To ensure optimal performance, use the swabs supplied in the test kit. Alternatively, sterile nylon, foam or rayon nasal swabs can be used for specimen collection. Do not use calcium alginate swabs.

#### Nasal swab

- To collect a nasal swab, insert a sterile swab into the nostril that presents the most secretion under visual inspection. If no secretion is visible, insert the sterile swab into the nostril that is most congested.
- Gently push the swab until resistance is met at the level of the turbinate (less than 2,5 cm into the nostril).
- Rotate the swab a few times against nasal wall.
- Slowly withdraw the swab while continuing to rotate it.

**Note:** In patients whose nasal cavity is dry, wet the swab with a sterilised physiological saline solution (not supplied in the test kit) in advance and then collect a sample with it.

#### Nasopharyngeal swab

- Insert the swab carefully into the nostril that presents the most secretion under visual inspection. Keep the swab near the septum, the floor of the nose while gently pushing the swab into the posterior nasopharynx.
- Rotate the swab several times.

#### Nasal wash

- With the patient's head hyper-extended, add sterile, normal saline to a nostril using a syringe. Use the minimal amount of saline which the procedure allows, as excessive volume may dilute the antigen in the specimen.
- To collect nasal wash, place a clean, dry specimen container directly under the nose with a slight pressure on the upper

lip. Tilt the head forward allowing the fluid to run out of the nostril into the specimen container.

- Repeat the procedure for the other nostril and collect nasal wash into the same specimen container.

**Note:** Normal saline, a syringe and a specimen container are not supplied in the test kit.

#### Nasal aspirate

- Insert an aspirating tube with a trap up to the depth of the nasal cavity. Connect another tube to the aspiration device applying a negative pressure to it. Aspirate nasal fluid into the trap.

**Note:** The aspiration device is not supplied in the test kit.

For nasal washes/aspirates, the volume of 1-3 mL is recommended. If transport medium is used, minimal dilution of specimens (~1 mL) is recommended.

#### Sample Storage & Transport

Specimens should be tested as soon as possible after collection. If sample transport is required, the following transport media, which have been tested and shown no interference with the test performance, are recommended: brain-heart infusion broth, Hanks' balanced salt solution, M5 media, saline or phosphate buffer solution.

Alternatively, samples may be stored refrigerated at 2-8°C or stored at room temperature (15-30°C) in a clean, dry, closed container for up to 8 hours prior to testing. Nasal washes or aspirates may also be stored frozen (-70°C or colder) for up to one month.

## 9. Test Procedure

**Bring the tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.**


1. For each specimen, remove a test cassette from the sealed foil pouch and place it on a clean and level surface. Label the test cassette with the patient or control identification. For the best results, the assay should be performed within one hour.
2. Gently mix the buffer and add 8 drops to the extraction tube.

#### 3a. Nasal/nasopharyngeal swab


- a) Insert the swab with collected nasal/nasopharyngeal specimen into the tube. Swirl the swab to mix the reagents and compress the walls of the tube against the swab.
- b) Remove the swab while compressing the walls of the tube against the swab. Discard the swab in accordance with your biohazard waste disposal protocol.

#### 3b. Nasal wash/aspirate

- a) Vortex or thoroughly mix the specimen. Do not centrifuge it, as the removal of cellular material may adversely affect test sensitivity.
- b) Transfer 300 µL of the specimen into the extraction tube, using a transfer pipette (not provided in the test kit).

4. Switch on the nal von minden Colibri and make sure that the test-specific RFID card is at hand.


**Note:** Every lot has its own calibration data. These are saved on a RFID card delivered with every test package. For the correct test interpretation, the lot-specific RFID card must be read by the nal von minden Colibri. Old RFID cards should be discarded after using up the tests with the corresponding lot to avoid confusion between the cards.

5. Adjust a dropper-cap to the extraction tube. Add 2 drops (approximately 100 µl) of the sample to the sample well (S) of the test cassette by gently squeezing the tube.

**Avoid trapping air bubbles in the specimen well (S).**

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

6. Start the timer or use the integrated timer in the nal von minden Colibri.  
7. Place the nal von minden Colibri onto the test cassette, so that it fits exactly into the cavity.  
8. Interpret test results after 15 minutes. Do not interpret test results after more than 20 minutes.


## 10. Result Interpretation

The nal von minden Colibri shows the qualitative result following an analysis on the display.

Note that this is a qualitative test only and it cannot determine the concentration of the analyte in the specimen.

### Invalid

The control line 'C' fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. The nal von minden Colibri will also display this result.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

**After the results have been interpreted, used tests should be discarded immediately in accordance with local regulations for potentially infectious materials.**

## 11. Calibration

The calibration data of the NADAL® Influenza A/B Scan test is saved on the included RFID card. It is therefore not necessary for the user to calibrate the device.

## 12. Quality Control

An internal procedural control is included in the test cassette: A coloured line appearing in the control line region 'C' is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

*Good laboratory practice (GLP)* recommends the use of control materials to ensure proper test kit performance.

## 13. Limitations

- The NADAL® Influenza A/B Scan Test is for professional *in-vitro* diagnostic use only and should only be used for the qualitative detection of influenza virus type A and/or B. The colour intensity of the test line should not be used for quantitative or semi-quantitative evaluation.
- Additional testing is required to differentiate any specific influenza virus A and B subtypes or strains in consultation with state or local public health departments.
- Both viable and non-viable influenza A and B viruses are detectable, using the NADAL® Influenza A/B Scan Test.
- The performance characteristics of the NADAL® Influenza A/B Scan Test have not been established for use in monitoring antiviral treatment or as a cell culture identification method.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- Monoclonal antibodies may fail to detect or detect with less sensitivity influenza A virus antigens that have undergone minor amino acid changes in the target epitope region.
- Failure to follow the sections 'Test Procedure' and 'Result Interpretation' may adversely affect test performance and/or invalidate the test result.
- Individuals who received nasally administered influenza A vaccine may have positive test results for up to three days after vaccination.
- Results obtained with this assay, particularly in the case of weak test lines that are difficult to interpret, should be evaluated in conjunction with other clinical information available to the physician.
- The etiology of respiratory infection caused by microorganisms other than influenza A or B virus cannot be established with this test.
- Children tend to shed virus more abundantly and for longer periods of time than adults. Therefore, testing specimens from adults may often yield lower sensitivity than testing specimens from children.
- Positive and negative predictive values are highly dependent on prevalence. False negative test results are more likely during peak influenza activity when the prevalence of disease is high. False positive test results are more likely during periods of low influenza activity when prevalence is moderate to low.

- A high-dose hook effect may occur where the colour intensity of test lines decreases as the concentration of antigen increases. If a hook effect is suspected, the dilution of specimens may increase the colour intensity of test lines.

#### 14. Performance Characteristics

##### Clinical Studies

A total of 280 nasal swabs from patients with suspected influenza A viral infection were tested, using the NADAL® Influenza A/B Scan Test. From these samples 108 were confirmed to be positive and 172 were confirmed to be negative by cell culture. The results are presented in the following table:

##### Nasal swabs for influenza A

		Cell Culture +	Cell Culture –	Total
NADAL® Influenza A/B Scan Test	Influenza A +	96	9	105
	Influenza A –	12	163	175
	Total	108	172	280

Positive agreement with cell culture: 96/108=88.9%

Negative agreement with cell culture: 163/172=94.8%

Total agreement with cell culture: (96+163)/280=92.5%

A total of 280 nasal swabs from patients with suspected influenza B viral infection were tested, using the NADAL® Influenza A/B Scan Test. From these samples 89 were confirmed to be positive and 191 were confirmed to be negative by cell culture. The results are presented in the following table:

##### Nasal swabs for influenza B

		Cell Culture +	Cell Culture –	Total
NADAL® Influenza A/B Scan Test	Influenza B +	73	7	80
	Influenza B –	16	184	200
	Total	89	191	280

Positive agreement with cell culture: 73/89=82%

Negative agreement with cell culture: 184/191=96.3%

Total agreement with cell culture: (73+184)/280=91.8%

A total of 190 nasopharyngeal swabs from patients with suspected influenza A viral infection were tested, using the NADAL® Influenza A/B Scan Test. From these samples 78 were confirmed to be positive and 112 were confirmed to be negative by cell culture. The results are presented in the following table:

##### Nasopharyngeal swabs for influenza A

		Cell Culture +	Cell Culture –	Total
NADAL® Influenza A/B Scan Test	Influenza A +	65	9	74
	Influenza A –	13	103	116
	Total	78	112	190

Positive agreement with cell culture: 65/78=83.3%

Negative agreement with cell culture: 103/112=92%

Total agreement with cell culture: (65+103)/190=88.4%

A total of 190 nasopharyngeal swabs from patients with suspected influenza B viral infection were tested, using the NADAL® Influenza A/B Scan Test. From these samples 85 were confirmed to be positive and 105 were confirmed to be negative by cell culture. The results are presented in the following table:

##### Nasopharyngeal swabs for influenza B

		Cell Culture +	Cell Culture –	Total
NADAL® Influenza A/B Scan Test	Influenza B +	70	6	76
	Influenza B –	15	99	114
	Total	85	105	190

Positive agreement with cell culture: 70/85=82.4%

Negative agreement with cell culture: 99/105=94.3%

Total agreement with cell culture: (70+99)/190=88.9%

A total of 300 nasal washes from patients with suspected influenza A viral infection were tested, using the NADAL® Influenza A/B Scan Test. From these samples 113 were confirmed to be positive and 187 were confirmed to be negative by cell culture. The results are presented in the following table:

##### Nasal washes for influenza A

		Cell Culture +	Cell Culture –	Total
NADAL® Influenza A/B Scan Test	Influenza A +	96	8	104
	Influenza A –	17	179	196
	Total	113	187	300

Positive agreement with cell culture: 96/113=85%

Negative agreement with cell culture: 179/187=95.7%

Total agreement with cell culture: (96+179)/300=91.7%

A total of 300 nasal washes from patients with suspected influenza B viral infection were tested, using the NADAL® Influenza A/B Scan Test. From these samples 94 were confirmed to be positive and 206 were confirmed to be negative by cell culture. The results are presented in the following table:

##### Nasal washes for influenza B

		Cell Culture +	Cell Culture –	Total
NADAL® Influenza A/B Scan Test	Influenza B +	82	9	91
	Influenza B –	12	197	209
	Total	94	206	300

Positive agreement with cell culture: 82/94=87.2%

Negative agreement with cell culture: 197/206=95.6%

Total agreement with cell culture: (82+197)/300=93%

A total of 180 nasal aspirates from patients with suspected influenza A viral infection were tested, using the NADAL® Influenza A/B Scan Test. From these samples 71 were confirmed to be positive and 109 were confirmed to be negative by cell culture. The results are presented in the following table:

**Nasal aspirates for influenza A**

		Cell Culture +	Cell Culture –	Total
NADAL® Influenza A/B Scan Test	Influenza A +	59	5	64
	Influenza A –	12	104	156
	Total	71	109	180

Positive agreement with cell culture: 59/71=83.1%

Negative agreement with cell culture: 104/109=95.4%

Total agreement with cell culture: (59+104)/180=90.6%

A total of 180 nasal aspirates from patients with suspected influenza B viral infection were tested, using the NADAL® Influenza A/B Scan Test. From these samples 48 were confirmed to be positive and 132 were confirmed to be negative by cell culture. The results are presented in the following table:

**Nasal aspirates for influenza B**

		Cell Culture +	Cell Culture –	Total
NADAL® Influenza A/B Scan Test	Influenza B +	41	5	46
	Influenza B –	7	127	134
	Total	48	132	180

Positive agreement with cell culture: 41/48=85.4%

Negative agreement with cell culture: 127/132=96.2%

Total agreement with cell culture: (41+127)/180=93.3%

**Analytical sensitivity/LOD**

The limit of detection (LOD) was determined by evaluating different concentrations of one strain of influenza A virus and one strain of influenza B virus, using the NADAL® Influenza A/B Scan Test. Multiple operators tested each concentration of the two influenza strains multiple times. The concentrations determined as the LOD levels for each strain tested are listed below.

Influenza A: A2/Aichi/2/68(H3N2),  $2.3 \times 10^3$  CEID<sub>50</sub>/test\*

Influenza B: Hong Kong 5/72,  $3.5 \times 10^3$  CEID<sub>50</sub>/test\*

\*CEID<sub>50</sub>: Chicken Embryo Infectious Dose

**Analytical reactivity**

The human, avian or other animal-derived influenza A and B virus strains listed below have been tested and showed a positive reaction with the NADAL® Influenza A/B Scan Test. Although specific influenza strains causing infection in humans may vary from year to year, they all contain the conserved nucleoproteins targeted by the NADAL® Influenza A/B Scan Test.

**Influenza virus strain**

A/Narita/1/2009 (H1N1)	A/Chicken/Yamaguchi/7/04 (H5N1)
A/NWS/33 10 (H1N1)	A/Chicken/Italy/99 (H7N1)
A/Hong Kong/8/68 (H3N2)	A/Chicken/Netherlands/03 (H7N7)
A2/Aichi/2/68(H3N2)	A/Swine/Hokkaido/2/81 (H1N1)
A/WS/33 (H1N1)	A/Duck/Tottori/723/80 (H1N1)
A/New Jersey/8/76 (HswN1)	A/Duck/Hokkaido/17/01 (H2N3)
A/Mal/302/54 (H1N1)	A/Duck/Mongolia/4/03 (H3N8)

A/Anhui/1/2013 (H7N9)	A/Duck/Czech Republic/56 (H4N6)
A/Shanghai/1/2013 (H7N9)	A/Duck/Pennsylvania/10128/84 (H5N2)
A/Hong Kong/156/97 (H5N1)	A/Turkey/Massachusetts/3740/65 (H6N2)
A/Hong Kong/483/97(H5N1)	A/Seal/Massachusetts/1/80 (H7N7)
A/Duck/Mongolia/119/2008 (H7N9)	B/Hong Kong 5/72
A/Duck/Mongolia/128/2008 (H7N9)	B/R5
A/Duck/Mongolia/147/2008 (H7N9)	B/Russia/69
A/Duck/Mongolia/129/2008 (H7N9)	B/Lee/40
B/Victoria/2/87 lineage	B/Yamagata/16/88 lineage
B/Maryland/1/59	

**Analytical specificity (cross-reactivity)**

To determine the analytical specificity of the NADAL® Influenza A/B Scan Test, 69 commensal or pathogenic microorganisms (24 viruses, 45 bacteria) that may be present in the upper respiratory tract have been tested. Positive and negative specimens were spiked with these microorganisms. Bacterial or yeast isolates were evaluated at a concentration of  $10^7$ - $10^8$  org/mL. Viral isolates were inoculated at a concentration of  $10^4$ - $10^8$  TCID<sub>50</sub>/mL. Adenovirus 18 and parainfluenza virus 3 were tested at  $10^2$  TCID<sub>50</sub>/mL. None of the microorganisms tested yielded a positive result with the influenza-negative samples or interfered with the detection of the influenza A or B positive samples. Both the negative and positive respiratory specimens were positive when spiked with influenza A strain A2/Aichi/2/68 (H3N2) or influenza B strain Hong Kong 5/72.

**Virus other than influenza A/B viruses**

Human adenovirus B, C	Adenovirus type 10, 18
Coxsackievirus A9, B5	Human herpesvirus 2, 5
Herpes simplex virus 1	Human rhinovirus 2, 14, 16
Mumps virus	Sendai virus
Respiratory syncytial virus	Rubella virus
Human coronavirus OC43	Echovirus 2, 3, 6
Measles virus	Parainfluenza virus 2, 3
Varicella-Zoster virus	

**Bacteria**

<i>Acinetobacter calcoaceticus</i>	<i>Moraxella catarrhalis</i>
<i>Bacillus cereus</i>	<i>Mycobacterium avium</i>
<i>Bacillus subtilis</i>	<i>Mycobacterium intracellulare</i>
<i>Bacteroides fragilis</i>	<i>Mycobacterium tuberculosis</i>
<i>Bordetella parapertussis</i>	<i>Mycoplasma pneumoniae</i>
<i>Bordetella pertussis</i>	<i>Neisseria meningitidis</i>
<i>Branhamella catarrhalis</i>	<i>Neisseria sicca</i>
<i>Chlamydia pneumoniae</i>	<i>Neisseria subflava</i>
<i>Citrobacter freundii</i>	<i>Nocardia asteroides</i>
<i>Corynebacterium diphtheriae</i>	<i>Proteus vulgaris</i>
<i>Enterobacter cloacae</i>	<i>Pseudomonas aeruginosa</i>
<i>Enterococcus faecalis</i>	<i>Serratia liquefaciens</i>
<i>Escherichia coli</i>	<i>Staphylococcus aureus</i>
<i>Gardnerella vaginalis</i>	<i>Staphylococcus epidermidis</i>

<i>Haemophilus influenzae</i>	<i>Streptococcus</i> Gruppe A, B, C, F, G
<i>Klebsiella oxytoca</i>	<i>Streptococcus mutans</i>
<i>Klebsiella pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Lactobacillus casei</i>	<i>Streptococcus salivarius</i>
<i>Lactobacillus plantarum</i>	<i>Streptococcus sanguis</i>
<i>Legionella pneumophila</i>	<i>Yersinia enterocolitica</i>
<i>Listeria monocytogenes</i>	

### Interfering substances

The following substances, normally present in respiratory specimens or that may be artificially introduced into the nasal or nasopharynx cavity, have been evaluated at the concentrations listed below. None of them were found to affect test performance of the NADAL® Influenza A/B Scan Test.

Substance	Concentration	Substance	Concentration
3 OTC nasal sprays	10%	Gualacol glyceryl ether	20 mg/mL
3 OTC mouthwashes	10%	Mucin	1%
3 OTC throat drops	10%	Mupirocin	250 µg/mL
4-acetamidophenol	10 mg/mL	Oxymetazoline	10 mg/mL
Acetylsalicylic acid	20 mg/mL	Phenylephrine	10 mg/mL
Albuterol	20 mg/mL	Phenylpropanolamine	20 mg/mL
Chlorpheniramine	5 mg/mL	Relenza® (zanamivir)	20 mg/mL
Dexamethasone	5 mg/mL	Rimantadine	500 ng/mL
Dextromethorphan	10 mg/mL	Tamiflu® (oseltamivir)	100 mg/mL
Diphenhydramine	5 mg/mL	Tobramycin	40 mg/mL
Doxylamine succinate	1 mg/mL	Triamcinolone	14 mg/mL
Flunisolide	3 mg/mL		

### Reproducibility study

A blind study, using the NADAL® Influenza A/B Scan Test was conducted at three separate clinical sites. Panels of blind-coded specimens containing negative, low positive (at the LOD) and moderate positive (above the LOD) influenza A and B viral samples were tested by three non-professional users per site to evaluate the reproducibility of the NADAL® Influenza A/B Scan Test. The users tested each sample multiple times on three different days. 96% of the samples tested showed the expected result.

### 15. References

1. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
2. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
3. CDC website: <http://www.cdc.gov/flu/>
4. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

Rev. 4, 2021-01-20 KBe/SDe

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consultense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico-diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límites de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

## Our Teams

### Germany:

#### Regensburg

Tel: +49 941 290 10-0  
Fax: +49 941 290 10-50

#### Moers

Tel: +49 2841 99820-0  
Fax: +49 2841 99820-1

#### Austria:

Tel: +49 941 290 10-29  
Free Tel: 0800 291 565  
Fax: +49 290 10-50  
Free Fax: 0800 298 197

#### UK & Ireland:

Tel: +49 941 290 10-18  
Free Tel – UK: 0808 234 1237  
Free Tel – IRE: 1800 555 080  
Fax: +49 290 10-50

#### France:

France Tel: 0800 915 240  
France Fax: 0800 909 493

#### Switzerland

Swiss Tel: 0800 564 720  
Swiss Fax: 0800 837 476

#### Belgium

Belgium Tel: 0800 718 82  
Belgium Fax: 0800 747 07

### Luxembourg

Lux. Tel: 800 211 16  
Lux. Fax: 800 261 79

### Spain:

Tel: +49 941 290 10-759  
Free Tel: 900 938 315  
Fax: +49 941 290 10-50  
Free Fax: 900 984 992

### Italy:

Tel: +49 941 290 10-34  
Fax: +49 941 290 10-50

### Poland:

Tel: +49 941 290 10-44  
Free Tel: 00 800 491 15 95  
Fax: +49 941 290 10-50  
Free Fax: 00 800 491 15 94

### Portugal:

Tel: +49 941 290 10-735  
Tel. Verde: 800 849 230  
Fax: +49 941 290 10-50  
Fax Verde: 800 849 229

### Netherlands:

Tel: +31 30 75 600  
Free Tel: 0800 0222 890  
Fax: +31 70 30 30 775  
Free Fax: 0800 024 9519

### Nordic countries:

#### Denmark

Tel: +31 703075 605  
Free Tel: 808 887 53

#### Finland

Tel: +31 703075 606  
Free Tel: 0800 918 263  
Free Fax: 0800 918 262

#### Norway

Tel: +31 703075 605  
Free Tel: 800 16 731

#### Sweden

Tel: +31 703075 605  
Free Tel: 020 79 09 06

#### Laboratory Diagnostics Team:

Tel: +49 941 290 10-40  
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Tel: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH