



Acidic softener for the thermal reprocessing of bedpans and urine bottles

thermosept® SKS

Our Plus

- prevents limescale
- foam-free, free of surfactants
- eco-friendly
- high material compatibility with stainless steel and different types of plastics

Application areas

thermosept® SKS is suitable for the thermal reprocessing of bedpans and urine bottles in any commercially available equipment. Product is applicable for all levels of water hardness.

Instructions for use

Dosage:

0.01 % per 1° dH (corresponds to 0.18 mmol/l water hardness), generally that corresponds to a dose of 0.1 - 0.3 % (1 - 3 ml/l water).

Dosing is performed by machine-integrated dosing pumps. Please check dosing pumps for acid resistance prior to application. If there is a large amount of limescale on the heating element descaling is recommended before using thermosept® SKS. thermosept® SKS prevents limescale in the steam generator, also with hard water. To avoid limescale in general, not only should regular maintenance be carried out but the canister should be replaced at regular intervals.



Product data

Composition:

hardness stabiliser, corrosion inhibitors, organic acids.

Chemical-physical data

Color	colourless
Density	ca. 1,08 g/cm ³ / 20 °C
Flash point	> 100 °C
Form	liquid
pH	ca. 1,3 / 20 °C
Viscosity, dynamic	No data available

Information for order

Item	Delivery form	Item no.
thermosept SKS SOKA 5 kg KA	1/Canister	on request
thermosept SKS 5 kg KA	1/Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Application aids

Application aids	Item no.
Can key for 5 + 10 l	135810

Environmental information

Schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: <http://www.schuelke.com/>.
For individual questions: Customer Sales Service Phone: +49 40 52100-666 E-Mail: info@schuelke.com



GEPRÜFTES
UMWELTMANAGEMENT
REG. NR. 0-150-00003



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

Schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com